
The Role of the CXCL16-CXCR6 Pathway in the Migration of Human-Induced **Spheroidal Neural Stem Cells in Vitro**



¹Department of Psychology and Neuroscience, University of North Carolina at Chapel Hill ²Division of Pharmacoengineering and Molecular Pharmaceutics, Eshelman School of Pharmacy, University of North Carolina at Chapel Hill


Introduction
 Glioblastoma (GBM) is an aggressive and incurable cancer of the central nervous system¹.
Median survival time: 15-23 months ² rate: 90% of patients ³
 There are numerous obstacles in treating GBM that have impeded the progress of developing effective therapies: X Rapid and extensive infiltration of surrounding brain tissue⁴ X Development of treatment resistance⁵ X Heterogenous tumor⁴ + microenvironment⁵ X Immunosuppressive microenvironment⁴ X Upregulation of efflux pumps⁴
 Human-induced spheroidal neural stem cells (hiNeuroSs) are a promising cell-based therapy for GBM due to their robust innate tumor-homing characteristics⁶. ✓ Can be engineered to continuously secrete biotherapeutics ✓ Allows for precise and persistent delivery of anti-tumor therapies to distant GBM tumor foci within 72 hours → reduces likelihood of tumor recurrence
Research Aim: Determine the mechanism behind hiNeuroS tumor-homing migration. Hypothesis: The CXCL16-CXCR6 chemokine pathway plays a role in the tumor-homing capabilities of hiNeuroSs.
 <u>CXCL16:</u> Chemoattractant Overexpressed in GBM^{7, 8} <u>CXCR6:</u> Receptor implicated in migration of GPCs⁹ and MSCs¹⁰
CXCL16 GBM CXCR6 hiNeuroS hiNeuroS hiNeuroS hiNeuroS Migration

Clara Mellows¹ & Dr. Shawn Hingtgen²

Methods

Agarose Gel Migration Assay:

1) 6-well plates prepared with 5mL agarose gel w/ ReNcell media per well. HiNeuroS-mCherry-FLuc seeded in channels on left side. On right side, 4 injections of either PBS or sCXCL16.



2) Fluorescent imaging using an EVOS M7000 imaging system @ 0, 6, 24, and 48-hours after treatment. 3 plates imaged with 3 images per well, per timepoint.



3) ImageJ software used to quantify hiNeuroS-mCherry signal in the region of interest, which started at the right-most signal in the 0-hour image for each location in the




plate. 4) Statistical analysis using Jamovi software: Repeated measures ANOVA with post-hoc Tukey.





Figure 1. Viability Assay for hiNeuroSs in Agarose Gel. Fluorescent hiNeuroS-mCherry signal intensity measured at 0, 24, and 48 hours after seeding 5x10⁴ or 1x10⁵ cells in agarose gel channels in the absence of any treatment. N = 3 images were quantified for each treatment group at each timepoint. * = p < 0.05, ns (not significant) = p > 0.05.





Discussion

- > Novel agarose migration assay developed to fit the unique characteristics of hiNeuroSs.
- > CXCL16 at higher doses (25 ng/mL) induced migration of hiNeuroSs over 48 hours *in vitro*.
- The CXCL16-CXCR6 axis may play a role in the tumor-homing capabilities of hiNeuroSs

Implications:

- Could inform further hiNeuroS optimization for the treatment of GBM
- Novel agarose gel migration assay; may be applicable for the study of other migrating cell lines if the appropriate media is substituted

Limitations:

- Variability in shape of channels and distance from treatment
- Possible cell death from lower cell count

Future Directions:

- Improve current study (increase cell count; increase time-course beyond 48 hours)
- Antagonize CXCR6 in hiNeuroSs \rightarrow Agarose migration assay
- CXCR6 constitutive expression in hiNeuroSs \rightarrow Agarose migration assay
- CXCR6 constitutive expression in hiNeuroSs \rightarrow *In vivo* migration studies

References

▣;;;;;□

Full Text PDF

Ask me about a full text PDF!

clmellows@hotmail.com

Acknowledgements

This project was supported by a David Bray Peele Memorial Research Award from the Department of Psychology and Neuroscience, University of North Carolina at Chapel Hill.